



MANITOBA RENAL PROGRAM

SUBJECT <ul style="list-style-type: none"> Setting up NxStage Cyclor for Home Hemodialysis treatment. 	SECTION 50.20 Home Hemodialysis
	CODE 50.20.08
AUTHORIZATION <ul style="list-style-type: none"> Professional Advisory Committee, Manitoba Renal Program 	EFFECTIVE DATE January 2018
	REVISION DATE

PURPOSE:

- To ensure the proper setup of the NxStage cyclor for the use of home hemodialysis.

POLICY:

- Registered Nurses in the Home Hemodialysis program who have received instruction and have demonstrated competency to the renal educator or delegate may instruct patient's on the Home Hemodialysis program on the proper technique in setting up the Nxstage cyclor for the propose of Home Hemodialysis.

EQUIPMENT:

- NxStage System One cyclor
- Two 10 mL sterile empty syringes
- Two 1 litre Normal Saline 0.9% sterile bags
- NxStage Cyclor Cartridge (CAR-172-C)
- NxStage PureFlow SL® and cabinet
- Medical tape
- Dry towel
- Wet towel
- Glass lens cleaner cloth

PROCEDURE:

- Perform hand hygiene.
- Gather and inspect supplies.
- Open the cyclor door.
- Turn on the power switch at the back of the cyclor.
- Wait until the yellow caution bars flash on the cyclor display area.

KEY POINT:

- Verify expiry dates, sterility and ensure accurate priming solution.
- Lift the cyclor door until it clicks and pull the handle to open the door fully. Ensure the handle remains up.
- Dust the blood leak detector sensor weekly using a glass-lens cleaning cloth.

PROCEDURE:

6. Remove the cartridge from the packaging. Make sure the happy face ☺ is facing up.
7. Tape the packaging/bag to the end of the table.
8. Remove the female connector (white plastic removable piece) from the cartridge and tape it to the top of the cyclor.
9. Remove the 2 connector tabs and seven pieces of tape.
10. Secure the cap connections – there are seven:
 - a. Arterial cap
 - b. Venous cap
 - c. Saline T line
 - d. Dialyzer port
 - e. Heparin port
 - f. Green cap
 - g. Yellow cap
11. Grab and pull entire tubing loop to the right.
12. Loosen tape and remove line organizer.
13. Ensure the happy face ☺ is facing up.
14. Take the cartridge and slide it into the cyclor from the top, resting the cartridge on the two metal bars inside the cyclor.
15. Hook the tubing on the cyclor IV pole.
16. Secure tubing into the cyclor air detectors.

There are 3 air detectors

 - a. Venous air detector –top left.
 - b. Arterial air detector –lower right.
 - c. Dialysate air detector – middle right.
17. Close the cyclor door.
18. Check the fluid sensor with the wet towel. Sensor should alarm. Dry immediately.
19. Remove one 1 litre saline bag from packaging. Keep the saline bag lying flat to spike.
20. Spike the saline bag with cartridge spike connector; use a rotating-twisting motion right to left.

KEY POINT:

- Use the packaging/ bag for all waste pre, during, post treatment.
- Do not remove the clear plastic tape on the dialyzer.
- Remember the air detectors are square.
- Before closing the cyclor door ensure all tubing is clear on both sides of the cyclor.
- Do not completely rotate in full circles this will tangle the lines.
- Do not throw cap of spike away, tape it to the top of the Cyclor. It is required later for disposing of the cartridge.

PROCEDURE:

21. Hang the saline bag on the cyclor IV pole.
22. Insert access pressure pod line into the access pressure pod port.
23. Press **"ADD FLUID"** button, this will initiate 23 minutes of cyclor priming.
24. After 23 minutes of priming, the cyclor will chime; the cyclor will display **"8's"** in the status bar and will provide an audible alarm.

Press **"MUTE"**

- The machine will show **"1234 AbCd eFgH"** in the status bar
- The machine will show **"987 654 321"** in the volume rate windows.

25. Press **"MUTE"**

26. The status bar will display:
"23.0 and 000"

27. Attempt to "flick (use the posterior aspect of the fingers)" to remove excess air bubbles within the system.
 - Start at the **"arterial red line"** (Starting at the priming spike at the saline bag) to the access pressure pod. Make sure to flick the access pressure pod too.
 - Remove the dialyzer from the cartridge tray by pulling and removing the plastic ribbon tape.
 - Flip the dialyzer over (blue clamp up) and gently rotate the dialyzer until there is a decrease in air bubbles.
 - Put the dialyzer in its holder with the **"blue clamp"** up. Make sure lines are not kinked.
 - Prime the post dialyzer port, then clamp **blue clamp** on the post dialyzer port.
 - Prime and clamp the heparin line as close to the blood lines as possible.
 - Flick (use the posterior aspect of the fingers) the venous patient line starting from the cartridge (top right corner) – this will chase the air back towards the saline bag.
 - Flick the saline line.

28. Press **"STOP"**.

KEY POINT:

- Turn entire slip tip of the access pressure pod line to the left until tight then secure with the locking collar by turning to the right.
- Perform patient assessment, vital signs, and calculate weight loss.
- The code verifies delivery system visually and audibly working.
- Confirms that Nxstage software is operating correctly.
- Do not press **STOP**.
- If **STOP** is pressed, the machine will require a re-prime for an additional 23 minutes.
- Flicking of dialysis lines ensures air removal.
- Setup may have to be flicked 2 -3 times to ensure air removal.
- Flick:



PROCEDURE:

29. Prime and close the white clamp on the saline (T) on the red arterial line
30. Drape the saline (T) onto the cyclor IV pole.
31. Find the thin line with green clamp and white clamp attached to the saline line spike.
32. Close the white clamp and the small green clamp on the line.
33. Loosen the cap from the SAK® chicken toe (keep sterile).
34. Disconnect the white and green lines from each other.
35. Attach the green line to green chicken toe line.
36. Attach the saline line to the saline (T) – White to White
37. Ensure that the green clamps (3) on the dialysate line are unclamped.
38. Ensure that the saline lines are clamped (small white clamps)
39. Close the two yellow clamps on the waste line – hanging on the cyclor IV pole.
40. Disconnect the yellow clamped waste line from the saline bag priming spike and connect to the yellow waste line on the PureFlow SL®.
41. Open the two yellow clamps.
42. Enter dialysate rate/hour (cyclor green bar display)
43. Enter UF rate/hour (cyclor yellow bar display)
44. Press the cyclor Volume toggle button. This will change the screen to the VOLUME screen. The cyclor red bar display will say "VOL".
45. Enter the amount of dialysate that will be used for the entire treatment in the cyclor green bar display.
46. Enter the amount of fluid removal expected over the course of the dialysis treatment in the cyclor yellow bar display.
47. Press the cyclor volume toggle button to return to the main screen or wait for approximately one minute and system will revert to the main screen.

KEY POINT:

- Keep lines sterile.
- Keep connections sterile.
- Keep connections sterile.
- Dialysate rate is a maximum 12 litres/hour. Total SAK® volume is 60 litres. Patient will need to change dialysate flow rate if dialyzing longer than 5 hours.

PROCEDURE:

48. Establish vascular access as per MRP policy #:
30.20.01 *Venipuncture of Arteriovenous Fistula/Graft* or
30.20.02 *Accessing and Locking Dialysis Central Venous Catheter (Anticoagulant/Thrombolytic/Antibiotic Locking)* or
30.20.03 *Establishing and Maintaining Buttonhole Venipuncture Sites.*
30.20.04 *Use of Closed Needleless Access Device with Hemodialysis Central Venous Catheters (CVC).*

49. Clamp arterial and venous blood lines on priming spike.

50. Perform hand hygiene. Change the half-filled normal saline bag to a full 1 litre normal saline bag if site determines necessary.

KEY POINT:

- Do not drop. Maintain sterility on all line connections.

- Changeover of saline bag ensures sufficient normal saline for rinseback and flushes during treatment.

DOCUMENTATION:

- Document cartridge serial number on treatment sheet.

REFERENCES:

Nxstage PureFlow SL® User Guide. Software version 1.13, 1.14, and 1.15. NC5342 Rev. D 2015-09-09

RESOURCES:

Occupational and Environmental Safety & Health – WRHA (HSC Unit). October 2017 Ergonomic Assessment